**Never stop looking®**

# SecondLook® CAD Metrics Addendum

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2.2.1 *Meme Ucu Saptamaları*..... 80

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2.2.3 *Pektoral Kas Ana Hattı*..... 80

This Page Intentionally left Blank.

## ENGLISH

### 1 Overview of Manual

This addendum describes all of the SecondLook® CAD Metrics. The CAD Metrics are to be used as adjunctive information and are not to be used as diagnostic information. Their intent is to provide the end user with further information about any given detection and or case based assessments.

### 2 CAD Metrics

The software of the SecondLook system computes basic metrics for each region of interest detected thus enabling compatible review workstations to present this information to users in a variety of ways. The CAD software version information is plainly visible on the review workstation, typically near the text indicating the number of CAD marks per image.

CAD Metrics are provided for microcalcification clusters and densities where applicable. All metrics are defined in the DICOM standard and the information is contained in the DICOM Mammography CAD Structured Report. All of these parameters can be configured on or off. Implementation is review workstation dependent. Please refer to your review workstation DICOM conformance standard for a list of compatible features.

These metrics are broken down into two categories. The detection metrics are specific to the individual detection(s). System features are defined as functionality available to the workstation should they be able to utilize it.

#### 2.1 *Detection Metrics*

##### 2.1.1 Microcalcification and Density Area Measurement

This metric is the area of the region that CAD analyzed. This measurement should not be considered the true lesion area of a region with a MassMark or CalcMark; rather, it is simply the area of the region that was analyzed (as shown by the mass outline or the microcalcification point outlines).

##### 2.1.2 Number of Microcalcification Points

This metric is the number of locations that CAD analyzed as microcalcification points in the cluster. This is not the true number of microcalcifications in the region, as CAD can both miss true microcalcification points as well as incorrectly analyze non-calcium regions as microcalcification points.

##### 2.1.3 Microcalcification and Density Longest Axis Measurement

This metric is the diameter of the region that CAD analyzed, measured by finding the maximum distance between any two points on the mass outline or microcalcification point outlines. This measurement should not be considered the true lesion size of a region with a MassMark or

CalcMark; rather, it is simply the diameter of the region that was analyzed (as shown by the mass outline or the microcalcification point outlines).

#### 2.1.4 Distance from Nipple Measurement

This metric is the distance from the center of the nipple to the center of the CAD mark. CAD also provides the locations of the two points that were used to make the measurement.

#### 2.1.5 Distance from Chest Wall Measurement

This metric is the distance from the chest wall to the center of the CAD mark. CAD also provides the locations of the two points that were used to make the measurement.

#### 2.1.6 Mass, Microcalcification Cluster and Individual Microcalcification Point Outlines

This is the boundary of the region or the individual microcalcification that CAD analyzed.

#### 2.1.7 CAD iNSIGHT Score

This metric is a score from 1-100 that indicates how strongly the CAD considers that it should mark a region. Regions with high scores will be marked by CAD before regions with low scores. This is not an assessment of probability of malignancy. The score can differ from a radiologist's assessment; malignant regions can mistakenly have low scores, while non-suspicious regions can mistakenly have high scores. A radiologist should not change their assessment based on this score.

## 2.2 System Features

### 2.2.1 Nipple Detections

The location of a reference point used by CAD that is often on or near the location of the patient's nipple.

### 2.2.2 Skin Line Outline

The outline of the breast tissue that CAD examined for suspicious regions. Areas outside of the outlined region were not examined by CAD.

### 2.2.3 Pectoral Muscle Outline

The outline of the region that CAD considered the pectoral muscle.

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