



IHE Integration Statement			
Vendor	Product Name	Version	Date
iCAD, Inc.	PowerLook	11.1	8/1/2022
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Mammography Image (MAMMO)	Evidence Creator	None	
Evidence Document (ED) Integration Profile	Evidence Creator	None	
Post-Processing Workflow (PWF) Integration Profile	Evidence Creator	None	
<u>Internet address for vendor's IHE information:</u>			
http://www.icadmed.com/assets/dtb185-ihe-integration-statement-powerlook-11.pdf			
Links to Standards Conformance Statements for the Implementation			
HL7	Not applicable		
DICOM	http://www.icadmed.com/assets/powerlook-11-dicom-conformance-statement.pdf		
Links to general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	