



<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
iCAD, Inc.	PowerLook	10.0	6/02/2020
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Mammography Image (MAMMO)	Evidence Creator	None	
Evidence Document (ED) Integration Profile	Evidence Creator	None	
Post-Processing Workflow (PWF) Integration Profile	Evidence Creator	None	
<b><u>Internet address for vendor's IHE information:</u></b> <a href="http://www.icadmed.com/assets/dtb153-ihe-integration-statement-powerlook-10.pdf">http://www.icadmed.com/assets/dtb153-ihe-integration-statement-powerlook-10.pdf</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	Not applicable		
<b>DICOM</b>	<a href="http://www.icadmed.com/assets/powerlook-10-dicom-conformance-statement.pdf">http://www.icadmed.com/assets/powerlook-10-dicom-conformance-statement.pdf</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	